Instructions for Use/Package Insert

Ortho Organizers® Ancor Pro® Temporary Orthodontic Anchor

Description: The *Ortho Organizers Ancor Pro* Temporary Orthodontic Anchor is a temporary self-drilling and self-tapping anchor. The threaded portion of the device is 1,6 mm in diameter and is available in lengths of 6 mm, 8 mm, and 10 mm. The supragingival head of the device has an upper button and undercut that serves as an attachment point for elastic chain, closed loop springs, and other auxiliary orthodontic devices. The tapered conical section below the button incorporates a transverse circular opening which also serves as an attachment point, accommodating wires up to 0,022" in diameter. Four vertical flats below the conical section enable use of tools to turn the device into its fully seated position.

Material: The *Ortho Organizers Ancor Pro* Temporary Orthodontic Anchor is made of Ti-6Al-4V titanium alloy.

Indications: The *Ortho Organizers Ancor Pro* Temporary Orthodontic Anchor is intended to provide a fixed anchorage point for attachment of orthodontic appliances to facilitate the orthodontic movement of teeth. *Ancor Pro* Anchors can be loaded immediately after placement. They can remain in place for the course of orthodontic treatment (up to 24 months).

Contraindications: **Absolute Contraindications:** Patients presenting with one or more of the following are not candidates for the *Ancor Pro* Temporary Orthodontic Anchor.

- Active, intra-oral infection
- Inadequate bone quantity or quality
- History of metal sensitivities
- History of metabolic disorders
- Active periodontitis
- Radiation therapy regimen

Relative Contraindications: Patients presenting with one or both of the following may be at higher risk for complications and should undergo placement of the *Ortho Organizers Ancor Pro* Temporary Orthodontic Anchor only after careful risk/benefit consideration.

- Drug, alcohol, tobacco abuse
- Inability to maintain good oral hygiene

Potential Adverse Effects and Complications:

The following potential complications are common to all devices of this type.

- Poor bone formation, osteoporosis, osteolysis, osteomyelitis, inhibited revascularization, or infection can cause loosening, bending, cracking or fracture of the device or premature loss of fixation with the bone, leading to nonunion.
- Migration, bending, fracture or loosening of the temporary orthodontic anchor.
- Metal sensitivity or allergic reaction to a foreign body.
- Pain, discomfort or abnormal sensation, due to the presence of the device.
- Increased fibrous tissue response around the fracture site and/or the temporary orthodontic anchor.
- Necrosis of bone.
- Inadequate healing.

Apart from these adverse affects there are always possible complications of any surgical procedure such as, but not limited to, infection, nerve damage, and pain which may not be related to the temporary orthodontic anchor.

Warnings:

- Use of the *Ortho Organizers Ancor Pro* Temporary Orthodontic Anchor should not be undertaken without specific training, experience, and thorough familiarity with orthodontic principles and procedures.
- Improper placement of the *Ortho Organizers Ancor Pro* Temporary Orthodontic Anchor may result in failure of fixation or damage to surrounding vital structures (such as roots of teeth). Preoperative radiographic evaluation must be conducted to ensure adequacy of underlying bone quality and absence of vital structures in intended implantation site.
- Avoid creation of an over-expose pilot hole when using the drill which may result in excessive mobility of the anchor caused by the removal of the stabilizing cortical bone.
- Use of dissimilar metals can accelerate corrosion leading to device failure of sensitization to constituent metals. To minimize corrosion risk, use of nickel titanium or other non-stainless steel archwires is strongly recommended.
- The *Ortho Organizers Ancor Pro* Temporary Orthodontic Anchor must never be reused. Even in the absence of visible damage, a temporary orthodontic anchor of any type is unable to support equivalent loads if used more than once.

Packaging and Sterilization:

The *Ortho Organizers Ancor Pro* Temporary Orthodontic Anchor is supplied clean and non-sterile and is intended to be sterilized prior to use. Each individual anchor, when purchased individually, is packaged in an autoclavable peel pouch contained within an outer autoclavable peel pouch. Steam sterilization following AAMI or ISO standards and a validated protocol to achieve a sterility assurance level (SAL) of 10⁻⁶ or better. Alternatively, dry heat sterilization may be used after first removing the device from the autoclavable pouches. As with steam sterilization, dry heat sterilization should be conducted using a validated protocol in accordance with relevant AAMI or ISO standards to a sterility assurance level (SAL) of 10⁻⁶ or better. The manufacturer's recommended sterilization options are:

- Autoclave at 121°C/250°F, gravity, for a minimum of twenty (20) minutes followed by a 30-minute dry cycle, or
- Dry heat at 160°C/320°F for 2 hours.

Additional sterilization instructions for individual parts (not received with the surgical kit):

- Individually ordered anchors and tools are shipped in an autoclave pouch.
- Inspect the autoclave pouch for punctures prior to steam sterilization. Use a new autoclave pouch (bag), if the autoclave pouch shipped with the product is damaged.
- If using the autoclave method, allow the pouch to cool for 5 to 10 minutes after sterilization, and check that the color-coded indicator has changed:

Blue — Brown

The changing of the color-coded indicator only indicates that the pouch has been exposed to a sterilization cycle.

Preoperative Procedure:

1. Acquire radiographs of the intended anchor location to determine adequacy of bone volume, bone density, space, and length of the orthodontic anchor.
 2. Record and pour impressions of both upper and lower arches.
 3. Create surgical guide and embed the *Ancor Pro* X-Ray Placement Guide (ref 500-026) in the desired position.
- Surgical Procedure:**
1. Have patient rinse mouth with 0.12% chlorhexidine gluconate for 30-60 seconds.
 2. Apply topical anesthetic to anchor insertion site. A local anesthetic can be injected in to the anchor insertion site if preferred.
 3. Position surgical guide intraorally and obtain a periapical radiograph of the site to ensure proper location of the *Ancor Pro* X-Ray Placement Guide (ref 500-026). If necessary, reposition the *Ancor Pro* X-Ray Placement Guide to the intended position.
 4. Drilling Technique (optional):
 - a. In areas of dense cortical bone, it may be desirable to drill through the soft tissue and penetrate the cortical bone if

drilling is necessary, use *Ancor Pro* Drill provided in the kit (ref 500-020).

- b. **Note:** The drill should only be used to create a slight recess (~0.5mm) in the cortical bone and should not be used to its full depth or diameter.

5. The self-drilling, self-tapping anchor may be placed either manually or with a handpiece.

- a. Manual Method:
 - i. Insert anchor into either the *Ancor Pro* Short Driver Insert (ref 500-022) or *Ancor Pro* Long Driver Insert (ref 500-023).
 - ii. Insert the driver into the *Ancor Pro* Thumb Wheel (ref. 500-024) or *Ancor Pro* Screw Driver (ref 500-027).
 - iii. Rotate clockwise, using mild to moderate pressure, to begin inserting anchor into the bone to facilitate the self-drilling, self-tapping feature of the anchor. If the anchor does not begin to become stable after a few complete rotations, slightly increase the amount of pressure being applied.
 - iv. Seat the anchor by placing the long driver Insert (ref 500-023) or short driver insert (ref 500-022) in either the thumbwheel or the manual screwdriver for anterior region. In difficult to access posterior region, use the *Ancor Pro* Ratchet (ref 500-025).
 - v. If significant resistance is encountered, unscrew the anchor and increase the depth of penetration using the drill.
 - vi. Final seating is achieved when the top of the anchor is at the desired tissue level.

* *Ancor Pro* Ratchet (ref 500-025) is not included in kit and needs to be ordered separately. During insertion step, the arrow on the *Ancor Pro* Ratchet should face away from the patient to ensure clockwise rotation.

- b. Handpiece Method:
 - i. Insert the anchor into the *Ancor Pro* Latch Lock Driver (ref 500-021).
 - ii. Using mild to moderate pressure at very low speed (10-20 RPM), begin inserting the anchor. If significant resistance is encountered, unscrew the anchor and increase the depth of penetration using the drill (ref 500-020)
 - iii. Final seating is achieved when the top of the anchor is at the desired soft tissue level.

Removal

If soft tissue has grown over the collar of the anchor, anesthetize the surrounding soft tissues using either a topical or local anesthetic.

2. Removal

- a. Manual Method:
 - i. Insert anchor into either the *Ancor Pro* Short Driver Insert (ref 500-022) or *Ancor Pro* Long Driver Insert (ref 500-023).
 - ii. Insert the driver into the *Ancor Pro* Thumb Wheel (ref 500-024) or *Ancor Pro* Screw Driver (ref 500-027). If necessary, use the *Ancor Pro* Ratchet (ref 500-025).
 - iii. Begin unscrewing the anchor by turning the driver in a counter-clockwise direction.
- * *Ancor Pro* Ratchet (ref 500-025) is not included in kit and needs to be ordered separately. During insertion step, the arrow on the *Ancor Pro* Ratchet should face the patient to ensure a counter-clockwise rotation.
- b. Handpiece Method:
 - i. Insert the *Ancor Pro* Latch Lock Driver (ref 500-021) into slow speed handpiece.
 - ii. With handpiece set in reverse, at very low speed (10-20 RPM), begin unscrewing the anchor until removal is complete.

Caution: Federal Law (USA) restricts this device to sale by or on the order of an orthodontist, oral surgeon, dentist or physician.

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© 2010 Ortho Organizers, Inc. All rights reserved.
Ancor Pro is a registered trademark of Ortho Organizers, Inc.
Ortho Organizers, Inc.
Manufactured by Ortho Organizers, Inc.