Distributed by:

Ortho Organizers, Inc.  
1822 Aston Avenue  
Carlsbad, CA 92008  
USA  
Tel: 800.547.2000, +(1) 760 448 8600  
Fax: 800.888.7244, +(1) 760 448 8607  
USASales@OrthoOrganizers.com  
IntlOrders@OrthoOrganizers.com

Visit the *Ortho Organizers'* Web site at  
**OrthoOrganizers.com**

**mdi Europa GmbH**  
Langenhagener Str. 71  
30855 Langenhagen, Germany  
49 511 39 08 9530

OMA Member  
Orthodontic Manufacturers  
Association

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**Storage:**  
These materials are designed to be used at room temperature of approximately 21° to 24°C (70° to 75°F). Shelf life at room temperature is two years for unopened kits and one year for opened kits when caps are replaced on the syringes immediately after dispensing.  
Unopened kits can be refrigerated at 4°C (40°F) to extend the shelf life.  
If stored in a refrigerator, allow material to acclimate to room temperature prior to use.  
Do not store materials near any eugenol-containing products. Do not expose materials to elevated temperature or intense light.

**Curing:**  
Light cure for 50 seconds from the occlusal and gingival direction. Curing from both occlusal and gingival will strengthen the set of the material. Cure thoroughly with light tip close to the material.  
If using the autocuring method, the material will set completely in 8 minutes.

**Application:**  
Place the mixed cement into the band and set the band.  
Clean off excess cement.  
Replace caps on the syringes.  
Mix Catalyst and Base completely for 10 seconds.

**Mixing:**  
NOTE: There is a 3-minute working time from the beginning of the mix.  
Dispense equal amounts of Lightform Catalyst and Lightform Base onto the mixing pad. Protect the dispersed pastes from light. Do not place the material under the operator light since this will greatly decrease the working time.

**Band/Tooth Preparation:**  
Roughen the inside of the band with a fine diamond to improve the adhesion of the cement to the band.  
Prepare teeth with a non-fluoride, oil-free prophyl paste. Wash with an oil-free air/water spray.  
Isolate teeth with cotton rolls or rubber dam.  
Blow dry completely with an oil/moisture-free air source.

**Lightform Band Cement Kit**  
REF 470-130, 470-131, 470-132, and 470-133  
INSTRUCTIONS FOR USE

## Lightform Bandzement Kit REF 470-130, 470-131, 470-132, and 470-133 GEBRAUCHSANWEISUNG

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