# Endura® Instrument Sterilization

## INSTRUCTIONS FOR USE

**Warning:** This product is NOT sterilized. If this product is to be used intraorally, it must be sterilized by an approved method prior to use.

Ortho Organizers tests have shown that **DRY HEAT STERILIZATION IS THE BEST METHOD TO PREVENT CORROSION.**

### Dry Heat: Dentrionx DDS 5000

- Instrument must be thoroughly cleaned!** Clean 5 minutes in an Ultrasonic General Purpose Cleaner. Use distilled or filtered water for cleaning. Make sure that the cleaner has a rust inhibitor in its makeup.
- Dry joint and tip thoroughly with compressed air or towel ensuring joint and tip are free from moisture.
- Place instruments onto Dentrionx racks. Begin 20-minute cycle. Upon completion of cycle, remove instruments.
- Lubricate instruments (joint) using **Silicone Lubricant Spray (Catalog #600-123/600-124).**

### Cox Rapid Heat Transfer Sterilizer

- Instrument must be thoroughly cleaned!** Clean 5 minutes in an Ultrasonic General Purpose Cleaner. Use distilled or filtered water for cleaning. Make sure that the cleaner has a rust inhibitor in its makeup.
- Dry joint and tip thoroughly with compressed air or towel ensuring joint and tip are free from moisture.
- Place instruments into Cox sterilization chamber. Begin 6 minute cycle. Upon completion of cycle remove instruments and place to be cooled.
- Lubricate instruments (joint) using **Silicone Lubricant Spray (Catalog #600-123/600-124).**

**NOTE: The following sterilization methods can create a high risk of corrosion if the procedures are not followed correctly!**

### Autoclave: National Sterile Quick 9000D

- Instrument must be thoroughly cleaned!** Clean 5 minutes in an Ultrasonic General Purpose Cleaner. Use distilled or filtered water for cleaning. Make sure that the cleaner has a rust inhibitor in its makeup.
- Dry joint and tip thoroughly with compressed air or towel ensuring joint and tip are free from moisture.
- Place instruments flat on Autoclave tray with jaws open. Use distilled water in sterilizer.
- Twenty minute cycle in the Autoclave (260°F and 25 psi). When equipment has cooled sufficiently, remove the instrument(s). Use compressed air or towel to ensure that joint and tip are free from moisture.
- Lubricate instruments (joint) using **Silicone Lubricant Spray (Catalog #600-123/600-124).**

**NOTE: Instruments left in Autoclave may show signs of corrosion if not removed when cycle is complete.**

### Chemiclave: Harvey 4000/Chemiclave

- Instrument must be thoroughly cleaned!** Clean 5 minutes in an Ultrasonic General Purpose Cleaner. Use distilled or filtered water for cleaning. Make sure that the cleaner has a rust inhibitor in its makeup.
- Dry joint and tip thoroughly with compressed air or towel ensuring joint and tip are free from moisture.
- Place instrument flat onto Chemiclave tray. Place layer of paper toweling between instruments. Leave jaw open. Use distilled water in sterilizer.
- 20-minute cycle in the Chemiclave (250°F and 25 psi). When cycle is complete depressurize the equipment. Allow instruments to cool sufficiently. Use compressed air or towel to ensure that joint and tip are free from moisture.
- Lubricate instruments (joint) using **Silicone Lubricant Spray (Catalog #600-123/600-124).**

**NOTE: Instruments left in Chemiclave may show signs of corrosion if not removed when cycle is complete.**

### Cold Sterilization: 2% Glutaraldehyde

- Instrument must be thoroughly cleaned!** Clean 5 minutes in an Ultrasonic General Purpose Cleaner. Use distilled or filtered water for cleaning. Make sure that the cleaner has a rust inhibitor in its makeup.
- Dry joint and tip thoroughly with compressed air or towel ensuring joint and tip free from moisture.
- Immerse instrument in **sterilization liquid for at least 10 hours.**
- Remove instruments from solution. Rinse in sterile water. Use compressed air or towel to ensure that joint and tip are free from moisture.
- Lubricate instruments (joint) using **Silicone Lubricant Spray (Catalog #600-123/600-124).**

**NOTE: This method of sterilization creates a high risk of corrosion for inserted instruments and is not recommended. Do not use solution containing Benzyl Ammonium Chloride (BAC) and/or Hydrogen Peroxide (H<sub>2</sub>O<sub>2</sub>). These chemicals are known to deteriorate the tungsten carbide inserts. Use of these chemicals will void the warranty.**

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## ENDURA® INSTRUMENT WARRANTY POLICY

### Warranty Policy:

Ortho Organizers warrants Endura instruments to be free from any defect in material or workmanship, corrosion, or separation or fracture of inserts for a period of one (1) year from the date of purchase. Ortho Organizers warrants Endura Mathieu instruments to be free from any defect in material or workmanship, corrosion, or separation or fracture of inserts for a period of ninety (90) days from the date of purchase. Failure to follow all recommended use, maintenance, care and cleaning procedures will void warranty claims. Proof of recommended maintenance may be required where applicable.

### Care and Cleaning:

All recommended care, cleaning and sterilization procedures, as described in these Instructions for Use, must be fully adhered to. Failure to follow the recommended procedures will void the warranty. The use of rust inhibitor and silicone lube in the instruments' joints is required after each cleaning and sterilization.

### Recommended Maintenance:

Warranty does not include normal routine maintenance. All instruments should be reconditioned and re-sharpened every six (6) months. Ortho Organizers will re-sharpen or replace at no cost any cutting pliers which fail to properly cut their intended wire sizes within the first six (6) months of purchase.

Except as expressly provided above, there are no other warranties on Ortho Organizers' instruments, expressed or implied, with regards to description, quality, merchantability, or fitness for a particular purpose. The liability for any warranty claim shall not exceed the purchase price of the instrument that is the basis of such claim.

### Warranty Returns:

Instruments which are defective under the terms of the warranty policy, and have not been misused or tampered with, may be returned to Ortho Organizers for repair or replacement at Ortho Organizers' option. Contact Ortho Organizers and obtain a Return Material Authorization (RMA) number prior to returning the instruments by calling Customer Service at 1-800-547-2000. Once an RMA number is assigned, ship the product with this RMA number visible on the outer shipping package, along with the original invoice or packing slip, to Ortho Organizers, Inc., 1822 Aston Avenue, Carlsbad, CA 92008-7306. Please note that an RMA number is mandatory, and any product received by Ortho Organizers without an RMA number visible on the packaging will not be accepted.

### Routine Maintenance and Non-Warranty Repairs:

Do not send instruments to Ortho Organizers for routine maintenance, sharpening, or re-conditioning six (6) months after purchase or later. Normal maintenance, including re-sharpening and reconditioning, may be obtained at:

PRECISION PLIER SERVICE  
(800) 229-9060 or (909) 590-2085  
14175 Telephone Ave, Suite D  
Chino, CA 91710  
www.precisionplierservice.com

Alternative maintenance providers may also be used. Ortho Organizers does not warrant actions by any maintenance service providers, and assumes no responsibility for the results of any maintenance services performed.

Distributed by:  
Ortho Organizers, Inc.  
1822 Aston Avenue  
Carlsbad, CA 92008  
USA  
Tel: 800.547.2000, +(1) 760 448 8600  
Fax: 800.888.7244, +(1) 760 448 8607  
USASales@OrthoOrganizers.com  
IntlOrders@OrthoOrganizers.com

**mdi Europa GmbH**  
Langenhagener Str. 71  
30855 Langenhagen, Germany  
49 511 39 08 9530

OMA Member  
Orthodontic Manufacturers  
Association