Distributed by:

Ortho Organizers, Inc.
1822 Aston Avenue
Carlsbad, CA 92008
USA
Tel: 800.547.2000, +(1) 760 448 8600
Fax: 800.888.7244, +(1) 760 448 8607
USASales@OrthoOrganizers.com
IntlOrders@OrthoOrganizers.com

Visit the *Ortho Organizers'* Web site at
OrthoOrganizers.com

mdi Europa GmbH
Langenhagener Str. 71
30855 Langenhagen, Germany
49 511 39 08 9530

OMA Member
Orthodontic Manufacturers
Association

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Plastic Bracket Conditioner
REF 470-170
MODE D'EMPLOI

IMPORTANT: KEEP REFRIGERATED
This Plastic Appliance Conditioner creates a chemical bond between any plastic bracket or acrylic appliance and the bonding cement. It can also be used to condition a composite restoration prior to bonding, or to rebond a bracket impregnated with adhesive.

BONDING PLASTIC APPLIANCES
Use a brush or sponge pellet and apply a thin, uniform coat of Plastic Bracket Conditioner to the plastic surface to be bonded. Wait at least 60 seconds before applying bonding cement. (If more than 10 minutes elapses, re-apply conditioner.)

BONDING TO A COMPOSITE RESTORATION
Use a green stone to lightly rough-up the composite to remove the shine. Clean area with air syringe. Using a brush or sponge pellet and apply a thin, uniform coat of Plastic Bracket Conditioner and allow to dry. If bonding cement is used, apply sealant to the conditioned composite surface and to the cement of the bracket. Place bracket on tooth and position.

REBONDING A BRACKET
Use a bur to roughen the hardened adhesive on the bracket base. With a sponge pellet, apply a thin, uniform coat of Plastic Bracket Conditioner and allow to dry. Apply a thin coat of sealant to the conditioned bracket base. Next, apply sealant to the etched enamel, then apply paste to the bracket, place bracket on tooth and position.

Plastic Edgewise Brackets Part #802-201-204 and 802-211-214.
When bonding Plastic Edgewise Brackets, Plastic Bracket Conditioner, Part #470-170, must be used.

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