

	Lower 2nd Molar	Lower 2nd Molar	
<b>Upper 2nd Molar</b>	Lower 1st Molar	Lower 1st Molar	
	Lower 2nd Bi	Lower 2nd Bi	<b>Upper 2nd Molar</b>
<b>Upper 1st Molar</b>	Lower 1st Bi	Lower 1st Bi	<b>Upper 1st Molar</b>
<b>Upper 2nd Bi</b>	Lower Cuspid	Lower Cuspid	<b>Upper 2nd Bi</b>
<b>Upper 1st Bi</b>	Lower Anterior	Lower Anterior	<b>Upper 1st Bi</b>
<b>Upper Cuspid</b>	Lower Anterior	Lower Anterior	<b>Upper Cuspid</b>
<b>Upper Lateral</b>	<b>Upper Central</b>	<b>Upper Central</b>	<b>Upper Lateral</b>
<b>Patient's Right</b>			<b>Patient's Left</b>


1872 Arroyo Avenue, Carlsbad, CA 92008-7306 Phone: 760 448 8600  
Fax: 760 448 8601 Internet: www.OrthoOrganizers.com

**ILLUSION PLUS™ Ceramic BRACKET IDENTIFICATION**

**COLOR-CODED SITE GUIDES ON OCCLUSAL WITH PALMER NOTIFICATION; TAB POINTS TO DISTAL.**




UPPER RIGHT CENTRAL - Red  
UPPER LEFT CENTRAL - Green




LOWER RIGHT CENTRAL - Yellow  
LOWER LEFT CENTRAL - Blue




UPPER RIGHT LATERAL - Red  
UPPER LEFT LATERAL - Green




LOWER RIGHT LATERAL - Yellow  
LOWER LEFT LATERAL - Blue




UPPER RIGHT CUSPIDS - Red  
UPPER LEFT CUSPIDS - Green



LOWER RIGHT CUSPIDS - Yellow  
LOWER LEFT CUSPIDS - Blue




UPPER RIGHT BICUSPIDS - Red  
UPPER LEFT BICUSPIDS - Green



LOWER RIGHT BICUSPIDS - Yellow  
LOWER LEFT BICUSPIDS - Blue

\*Identical Brackets  
\*\*Identical Brackets



COLOR-CODED SITE GUIDES ON OCCLUSAL WITH PALMER NOTATION; TAB POINTS TO DISTAL.

**NOTE: These brackets are warranted for single use only.**  
Made in USA.

INS700  
082505