







ITEM	DESCRIPTION	CATALOG NO.																												
	<p>Extraoral Laser Welded Facebow (Pro-Bow™) This facebow is comfortable for the patient, plus it is easy to adjust and durable. The laser welded joint is much smoother and cleaner than a solder joint. The bow is made of solid wire with no hollow wire inserted. Hygienic, corrosion resistant, and more visually attractive. Permanently marked with a centerline. Wire diameter: .045; Angulation: 0°. Vertical loops on inner bow are adjustable.</p> <table border="1" data-bbox="427 646 1097 926"> <thead> <tr> <th>INNER ARCH</th> <th>SIZE NUMBER</th> <th>MM LENGTH</th> <th>CATALOG NO.</th> </tr> </thead> <tbody> <tr> <td>.045</td> <td>1</td> <td>83</td> <td>425-011</td> </tr> <tr> <td>.045</td> <td>2</td> <td>90</td> <td>425-012</td> </tr> <tr> <td>.045</td> <td>3</td> <td>97</td> <td>425-013</td> </tr> <tr> <td>.045</td> <td>4</td> <td>104</td> <td>425-014</td> </tr> <tr> <td>.045</td> <td>5</td> <td>111</td> <td>425-015</td> </tr> <tr> <td>.045</td> <td>6</td> <td>Straight</td> <td>425-016</td> </tr> </tbody> </table>	INNER ARCH	SIZE NUMBER	MM LENGTH	CATALOG NO.	.045	1	83	425-011	.045	2	90	425-012	.045	3	97	425-013	.045	4	104	425-014	.045	5	111	425-015	.045	6	Straight	425-016	
INNER ARCH	SIZE NUMBER	MM LENGTH	CATALOG NO.																											
.045	1	83	425-011																											
.045	2	90	425-012																											
.045	3	97	425-013																											
.045	4	104	425-014																											
.045	5	111	425-015																											
.045	6	Straight	425-016																											
	<p>Safety Neck Pad These modules greatly reduce the risk of slingshot rebound from headgear. Reassembly after release is easy. The stainless steel coil spring allows the doctor to apply specific forces depending on which hole he inserts the facebow into.</p> <p>Neck Pad & Safety Release with Stainless Steel Spring</p> <ul style="list-style-type: none"> • SS Black Safety Neck Pad (Black Clip) • SS Black Safety Neck Pad (Glow-in-the-Dark Clip) 	<p>426-002 426-003</p>																												
	<p>Safety Neck Pad Components</p> <ul style="list-style-type: none"> • Black Pad (5 per pack) • Black Clip (10 per pack) • Glow-in-the Dark Clip (10 per pack) 	<p>426-000 426-020 426-021</p>																												
	<p>Velpro™ Safety Neck Pad Soft <i>Velpro</i> Material offers maximum patient comfort, while greatly reducing slingshot rebound from headgear. Reassembly after release is easy. The stainless steel coil spring allows the doctor to apply specific forces depending on which hole he or she inserts the facebow into. Comes as a complete assembly only with black clip.</p> <p>Black Purple Pink Blue Teal</p>	<p>426-070 426-071 426-072 426-073 426-074</p>																												

ITEM	DESCRIPTION	CATALOG NO.
	<p>Face Mask (Reverse Pull) For use in Class III cases, made of stainless steel and durable plastic. Features a lightweight, comfortable, design. Plus, comes in a convenient carrying case. Easily adjusts to fit all patients.</p>  <p>Blue Red Teal Pearl</p>	<p>426-055 426-056 426-057 426-058</p>
	<p>Extraoral Elastics Available only in amber (1000 per pack) 8 oz 16 oz</p>	<p>408-073 408-074</p>
	<p>Safety High-Pull Head Caps Extra heavy duty construction for sturdy dependability. Easy to adjust.</p> <p>SS Black High Pull Cap</p> <ul style="list-style-type: none"> • Black Clip • Glow-in-the-Dark Clip <p>Safety High-Pull Head Cap Components</p> <ul style="list-style-type: none"> • Black Caps (5 per pack) • Black Clips (10 per pack) • Glow-in-the-Dark Clips (10 per pack) 	<p>426-012 426-013 426-010 426-020 426-021</p>
	<p>Class II Bow™ Designed for use with Class II elastics to distally drive or maintain anchorage of upper molars. Patient acceptance is much higher with <i>Class II Bow</i> over conventional headgear. Pre-soldered with pliable brass hooks for elastic use. Adjusting loops may be formed in wire, if desired. .045" wire diameter. (10 per pack)</p>	<p>100-090</p>
	<p>Chin Cup Designed for use with either high-pull or vertical-pull headgear. Felt lined for comfort. Rigid plastic cup will not squeeze patient's chin or cheeks.</p>	<p>426-050</p>