



Triple Play™

Clear Aligner Trays

Only Three Trays to a Picture Perfect Smile

Triple Play Clear Aligner Trays from the *O² Orthodontic Laboratory™* offer a comfortable, convenient and aesthetic solution to correct minor anterior crowding and rotations in the upper and lower arches. In addition, *Triple Play* is an ideal solution for post-retention relapse in cases that have already undergone more extensive orthodontic treatment. With a five day in lab turn-around, treatment can be started almost immediately and be completed in approximately six weeks.

Detailed Custom Analysis


Comprehensive analysis from our proprietary software system will provide patient specific case planning. The analysis will clearly show the number of trays needed and, unlike other solutions on the market, will note the exact amount and locations of any necessary interproximal reduction (IPR). Your case will receive detailed one-on-one attention from our experienced lab specialists, and support is only a call or e-mail away!

Cost Effective

With no set fee, you only pay for what each case requires – one, two, or a maximum of three trays depending on the severity of the case. Each value-priced aligner set-up includes model duplication, resetting the teeth and the actual tray.

Minimum Commitment

Only requiring a maximum of three sequential trays per arch, the time and cost commitment of *Triple Play* Clear Aligners is minimal. In addition, there are no contracts or special certifications required to utilize this system.



Sample of a comprehensive detailed custom analysis.

Auxiliary Products

Ortho Organizers® offers the following products to assist with interproximal reduction. These products give you the exact space you need and help shape the interproximal area.

Mini-Stripper™ Kit – easy to control in tight spaces and safer to use than other interproximal stripping devices. Single-sided strips prevent accidental stripping of adjacent tooth. Composed of Stainless Steel, all kit components are autoclavable up to 275°F/135°C. Catalog No. 4150-583



Separating Abrasive Strips – made of durable diamond abrasive, each strip has been electronically bonded to surgical grade steel base to resist stretching, breaking and corrosion. Strips are .003" thick and available in fine, medium or coarse in 4mm, 6mm or 8mm. (10 per pack)



Size	Fine	Medium	Coarse	Double-Sided Medium
4mm	106-220	106-221	106-222	106-221D
6mm	106-223	106-224	106-225	106-224D
8mm	—	106-227	106-228	106-227D

Abrasive Strip Holder – holds strips firmly to ensure better control during teeth shaping. Catalog No. 106-215



Call today, **888.851.0533**, or email **O2OrthoLab@OrthoOrganizers.com**, for more information on *Triple Play* Clear Aligner Trays. Don't forget to ask us about free inbound 2-day shipping via UPS and the laboratory starter kit.



For more information on our products, orthodontic laboratory services and educational offerings, please contact us:
 Customer Service: 888.851.0533
 In the U.S.: 800.547.2000
 To fax an order: 800.888.7244
 To email an order: USASales@OrthoOrganizers.com



A HENRY SCHEIN® COMPANY
 1822 Aston Avenue
 Carlsbad, CA 92008
 USA
OrthoOrganizers.com